

20th Annual Meeting of the TSE EURL
Rome, 2nd – 3rd October 2023

2022 round of TSE EURL EQAs:
Results feedback



**Confirmatory BSE and
Scrapie Immunoblotting**

Pier Luigi Acutis
IZSPLVA - Turin



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome

NRLs participating in Proficiency Testing 2022

23 and 24 laboratories for Confirmatory BSE and Scrapie



Austria, Belgium, Bulgaria, Croatia, Czech Republic, Finland, France, Germany, Greece, Hungary, Iceland*, Ireland, Italy, The Netherlands, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland and United Kingdom

** Only Confirmatory Scrapie*



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome



PLANNING

Each NRL has been identified with a new code, different for the two PTs

Each round included 5 samples of frozen brainstem homogenates of cattle and sheep

Samples were stored at -80 °C until shipment by the courier on dry ice

Set of samples included information on their preparation, manipulation and transmission results



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome



Shipping and Transmission results


Samples were shipped
on 21st and 22nd
November 2022



TSE EURL
Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome

**PROFICIENCY TESTING FOR VETERINARY LABORATORIES
GUIDELINES FOR FILLING THE RESULTS FORM**

Version num 2- Nov2020



Deadline for
submission of
results was on 16th
December 2022



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome



Transmission of the results

- All NRLs submitted the results of both PTs by the deadline.
- The submission included the raw data (pictures of the blots, readings), the information about the kit batches used for the analysis plus any information considered important for the correct interpretation of the results



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome



EURL EQA round 2022: Quality Control



Confirmatory BSE and Scrapie Immunoblotting

ISO 17034

Requirements for Producers of

Reference Materials



QUALITY



Homogeneity

Stability



ISO/IEC 17043
Conformity assessment – General requirements for proficiency testing



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome

BSE CONFIRMATORY IMMUNOBLOTTING

The samples provided for the TSE quality assessment are homogenates of brainstem material from a panel of BSE-positive and BSE-negative animals

The tissue have been homogenised using a 50/50 mix of tissue and distilled water

Sample No	Result	Interpretation
WB2201	negative	////
WB2202	positive	medium
WB2203	positive	strong
WB2204	negative	////
WB2205	negative	////

Sample WB2201 = sample WB2204 = sample WB2205



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome

SCRAPIE CONFIRMATORY IMMUNOBLOTTING

The samples provided for the TSE quality assessment are homogenates of brainstem material from a panel of Scrapie-positive and Scrapie-negative animals

The tissue have been homogenised using a 50/50 mix of tissue and distilled water

Sample No	Result	Interpretation	Characterisation
WS2201	positive	strong	classical scrapie
WS2202	positive	medium	classical scrapie
WS2203	positive	medium	classical scrapie
WS2204	negative	////	///
WS2205	positive	strong	classical scrapie

Sample WS2201 = sample WS2205; sample WS2202 = sample WS2203



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome

Processing of results

Elaboration of a report with the PTs results for
BSE and Scrapie round



Correct results



laboratories suitable to
perform confirmatory tests

Incorrect results



analysis on a new
set of samples



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome

Results Tabulation for BSE Confirmatory Immunoblot Scheme: WB22

Primary Tests	No. NRLs
TeSeE Western Blot - Bio-Rad	19
In house method	2
OIE SAF immunoblot	1
PrioWestern BSE Kit - ThermoFisher Scientific	1



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome

Results Tabulation for BSE Confirmatory Immunoblot Scheme: WB22



**TeSeE Western Blot Bio-Rad used four batches
1H0048 – 2C0049 – 2H0050 – 2F0046**

ThermoFisher
SCIENTIFIC

**PrioWestern BSE Kit – ThermoFisher Scientific used one batch
W200201G**



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome



BSE confirmatory Immunoblot – Scheme WB22




All laboratories have passed this PT round

Critical issues


TeSeE WB




TeSeE WB



TeSeE WB



PrioWestern



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome

Results Tabulation for Scrapie Confirmatory Immunoblot Scheme: WS22

Primary Tests	No. NRLs
TeSeE Western Blot Bio-Rad	16
In house method	3
APHA Bio-Rad Hybrid §	2
OIE SAF immunoblot	1
TeSeE Western Blot Bio-Rad *	1
PrioWestern BSE Kit - ThermoFisher Scientific**	1

§ APHA Bio-Rad Hybrid – used 2 mAbs: Sha31 and P4

* TeSeE Western Blot: detection also with P4

**PrioWestern BSE Kit - used 2mAbs: 6H4 and P4



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome

Results Tabulation for Scrapie Confirmatory Immunoblot Scheme: WS22



**TeSeE Western Blot Bio-Rad used four batches
1H0048 – 2C0049 – 2H0050 – 2F0046**

ThermoFisher
SCIENTIFIC

**PrioWestern BSE Kit – ThermoFisher Scientific used one batch
W200201G**



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome






Scrapie confirmatory Immunoblot – Scheme WS22



All laboratories have passed this PT round

Critical issues



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome





Scrapie confirmatory Immunoblot – Scheme WS22



All laboratories have passed this PT round

Critical issues



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome

Confirmatory BSE and Scrapie Immunoblotting

The forthcoming 2024 round of EQAs Confirmatory BSE and Scrapie Immunoblotting



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome

About shipping

- The shipment of the samples is scheduled in June
- BSE and Scrapie confirmatory PTs will be sent along with sample sets for the rapid test
- The samples will be packed with dry-ice and the shipment will be scheduled at the beginning of the week
- E-mail to communicate the shipment

Please transmit the import permit as soon as possible, if required and expired



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome

Submission PTs results

- The submission of results is expected by 3 weeks after the shipment
- The results submitted after the deadline will not be evaluated
- The results will be submitted by EURL web site
- Please include in the submission the raw data (pictures of the blots, readings) and all the information about the kit batches used for the analysis
- Report any information considered important for the correct interpretation of the results

**No extension of the
deadline will be made**



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome

PTs results

**The results of the PTs round 2023 and 2024
will be reported at the next TSE EURL Meeting**



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome

EURL Contact Immunolotting



Pier Luigi Acutis

e-mail address: pierluigi.acutis@izsto.it

telephone n.: +39 011 2686324

Maria Mazza

e-mail address: maria.mazza@izsto.it

telephone n.: +39 011 2686282



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome

THANKS FOR YOUR ATTENTION



TSE EURL

Istituto Zooprofilattico Sperimentale del Piemonte, Liguria e Valle d'Aosta – Turin
Istituto Superiore di Sanità - Rome