


# **EURL-AP Activities**

Joint Workshop EURL-TSE/AP - Rome, Italy, 14-15 May 2024



#### **Missions for EURLs**







### **Specificity of EURL-AP**




Light microscopy

Real-time PCR

Mass spectroscopy UHPLC-MS/MS









Challenge for the EURL-AP and the NRLs:
To become experts of very different techniques



### **Specificity of EURL-AP**







2<sup>nd</sup> challenge : sample and legislation complexity

adequate combination of techniques



relevant conclusion







### **Combined proficiency tests**



- Organisation of annual PT combining microscopic and PCR analyses
- Under ISO17043
- Participation of 3<sup>rd</sup> country official labs Iceland, Serbia, Norway China, Peru, Argentina, Australia, South Africa

Samples = challenge to produce Combination of methods not anymore under evaluation

| _aborato     | ory identificati                          | on code :   | 1                    |                           |   |            |                       | Cenn.                |                 |
|--------------|---|---|----------------------|---------------------------|---|------------|-----------------------|----------------------|-----------------|
|              |   | v agreement was   | not signed           |                           |   |            |                       | Antine               | al <b>Pro</b> l |
| Sampl<br>e N | CROSCOPY AN<br>Terrestrial<br>vert. part. | ALYSES: Details of ter pa   |                      | Terrestrial invert. part. | Details of terrestrial invert.<br>Part. | Fish part. | Details of fish part. | Sedimen<br>tation(s) |                 |
|              |   |   |                      |                           |   |            |                       |                      |                 |
|              |   |   |                      |                           |   |            |                       |                      |                 |
|              |   |   |                      |                           |   |            |                       |                      |                 |
|              |   |   |                      |                           |   |            |                       |                      |                 |
| Free cor     | nment                                     |   |                      |                           |   |            |                       |                      |                 |
| Free cor     |   |   |                      |                           |   |            |                       |                      |                 |
|              | LYSES:<br>Cut-off at 15 (.                | <i>5/orpig</i> ) oopies :   | Ruminant             | Poultry                   | Pigoyoles                               |            |                       |                      |                 |
|              | LYSES:<br>Cut-off at 15 (.<br>Copy numb   | <i>5/br.p(g</i> ) copies :<br>set at the out-off :<br>Master mix used : | Ruminant             | Poultry                   |   |            |                       |                      |                 |
|              | LYSES:<br>Cut-off at 15 (.<br>Copy numb   | er at the out-off :   | Ruminant Poultry DNA | Poultry                   | cycles                                  | Free comm  | ent                   | _                    |                 |
| PCR ANA      | LYSES:  Cut-off at 15 (.  Copy numb       | per at the cut-off:<br>Master mix used:<br><b>Pig</b>                   | Poultry              | Poultry                   | cycles                                  | Free comm  | ent                   | _                    |                 |
| PCR ANA      | LYSES:  Cut-off at 15 (.  Copy numb       | per at the cut-off:<br>Master mix used:<br><b>Pig</b>                   | Poultry              | Poultry                   | cycles                                  | Free comm  | ent                   |                      |                 |
| PCR ANA      | LYSES:  Cut-off at 15 (.  Copy numb       | per at the cut-off:<br>Master mix used:<br><b>Pig</b>                   | Poultry              | Poultry                   | cycles                                  | Free comm  | ent                   |                      |                 |



## **Working groups**



- Light microscopy
- PCR
- Mass spectroscopy
- Hydrolysed proteins
- GTH

### Expertise and support to:

- Choice of strategy to follow (e.g. PCR targets)
- Redaction of protocols (SOPs see after)





## **Standard Operating Procedure (SOPs)**



- Complements to the regulations
- Technical support,
- explanations to the Annex VI
- Better flexibility in case of changes in the protocols



European Union Reference Laboratory for Animal Proteins in feedingstuffs

**EURL-AP Standard Operating Procedure** 

#### Detection of ruminant DNA in feed using real-time PCR

Method development: TNO Triskelion bv Method assessment and validation: EURL-AP

| Version 1.0   | Last major revision |  |
|---|---------------------|--|
| Alessandro BENEDETTO Gilbert BERBEN Hermann BROLL Dilvier FUMIËRE Christoph HALDEMANN Lotte HOUGS Aline MARIEN Ingrid SCHOLTENS | EURL-AP team        |  |

Version 1.3

Page 1 on 16


Publication date 20.05.2021 Applicable on 25.05.2021



# **Standard Operating Procedure (SOPs)**



### https://www.eurl.craw.eu/





Missions

Our mission, in the global framework of the eradication of transmissible spongiform encephalopathies (TSE), focuses on the detection and identification of prohibited animal probability of the control o



Legal sources and SOPs

In Europe a very strict legislation has been set up over the years. This resulted in the "total feed ban" or the total ban of PAPs in feed. Each member state of the European Union needs



Publications

Linked to its activities, the EURL-AP team is publishing numerous papers and reports towards the scientific community or more accessible information dedicated to the general public.



Micrograph Collection

Access here to the EURL-AP micrograph collection, the largest on-line repository of reference pictures for animal remains identification, a tool supporting the feed ban control.



## Micrograph collection


























Picture of Black Soldier Fly (*Hermetia illucens*) larvae cuticule

- Started in 2007;
- Basically intended to help operators;
- Over 900 high quality records;
- Public tool intended for everyone and available on the website;
- Free use meanwhile a courtesy mention.



### **Development of E-tools**



#### On-line training



### Microscopy training

 Access to a virtual microscope that allows the participants to observe slides named according to their composition



#### Collaboration with the JRC



#### Production of calibrants for the PCR methods



ERM-AD482 RUMINANT pDNA CALIBRANT



ERM-AD483 PORCINE pDNA CALIBRANT



ERM-AD484k POULTRY pDNA CALIBRANT



## **Unexpected action**



#### Support to the EU Commission in the Horse meat crisis

- Not in the scope of the EURL-AP but
- EURL-AP sea mammalian PCR test not fully specific



→ development of a horse PCR test

- Horse PCR test at disposal of the control labs on request of the Commission
- Redaction of protocols
- Production of "reference" material for semi-quantitation



# **Acknowledgements to**





The EURL-AP team...



...and its NRLs network

