

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 1 di 15	DATA:



TSE EURL
Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta - Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RB24	Report No. 01/2025	Distribution date: 02 nd -03 rd December 2024
--	--------------------	---

Table 1. General information about 'routine' methods applied by the participants in the RB24 PT

NRL ID code	Test Method used	Kit Batch Number	Kit Expiry Date	Positive Control 1 result	Positive Control 2 result	Positive Control 3 result	Positive Control 4 result	Negative Control 1 result	Negative Control 2 result	Negative Control 3 result	Negative Control 4 result	Positive Cut Off value	Methods details
022	TeSeE SAP COMBI	4J0094	02/05/2025	1.721	1.688			0.006	0.006	0.008	0.007	0.217	
027	TeSeE SAP COMBI	4J0094	02/05/2025	2.492	2.393			0.009	0.010	0.009	0.007	0.219	
091	TeSeE SAP COMBI	4H0095	04/07/2025	3.249	3.268			0.007	0.011	0.008	0.008	Negative cut-off < 0.197; Positive cut-off >= 0.219	
132	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	BY445	28/02/2025	3.324	3.896			0.025	0.023	'-	'-	0.144	
135	IDEXX HerdChek BSE/Scrapie Antigen Test-ULTRASHORT PROTOCOL	BY445	28/02/2025	2.453	2.494			0.029	0.032	'-	'-	0.151	
181	IDEXX HerdChek BSE/Scrapie Antigen Test-ULTRASHORT PROTOCOL	EY700	05/06/2025	3.999	3.999			0.024	0.022	0.024	0.022	0.143	
192	IDEXX HerdChek BSE/Scrapie Antigen Test-ULTRASHORT PROTOCOL	AY185	08/01/2025	3.999	3.999			0.020	0.021	'-	'-	0.140	
291	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	AY185	08/01/2025	2.325	2.360			0.019	0.019			0.139	
326	IDEXX HerdChek BSE/Scrapie Antigen Test-ULTRASHORT PROTOCOL	EY700	05/06/2025	2.709	2.902			0.014	0.017	Na	Na	0.136	
329	TeSeE SAP COMBI	4C0093	14/01/2025	3.649	3.633			0.010	0.010	0.011	0.006	0.197	
358	TeSeE SAP COMBI	4J0094	02/05/2025	3.445	3.098			0.013	0.015	0.015	0.013	0.224	

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 2 di 15	DATA:



TSE EURL
Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta - Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RB24	Report No. 01/2025	Distribution date: 02 nd -03 rd December 2024
--	--------------------	---

390 433 (*)	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	CY816	21/03/2025	2.923 / 2.947	2.896 / 2.905			0.026 / 0.026	0.024 / 0.021	'-	'-	0.146 / 0.142	
051 397 (*)	IDEXX HerdChek BSE/Scrapie Antigen Test-ULTRASHORT PROTOCOL	EY700	05/06/2025	3.562	3.542			0.019	0.019			0.139	
415	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	EY700	05/06/2025	4.173	3.441			0.017	0.016	0.018	0.019	0.137	
595	TeSeE SAP COMBI	4C0093	14/01/2025	3.405	3.385			0.005	0.009	0.009	0.008	0.218	
596	TeSeE SAP COMBI	4C0093	14/01/2025	3.351	3.429			0.007	0.008	0.007	0.008	0.218	
651	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	AY185	08/01/2025	3.316	3.368			0.027	0.026	0.027	0.026	0.146	
672	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	EY700	05/06/2025	3.484	3.481			0.021	0.021	0.021	0.021	0.141	
703	IDEXX HerdChek BSE/Scrapie Antigen Test-ULTRASHORT PROTOCOL	CY816	21/03/2025	>4 (overflow)	>4 (overflow)			0.024	0.022	N/a	N/a	0.120	background corrected OD reported for all samples
740	IDEXX HerdChek BSE/Scrapie Antigen Test-ULTRASHORT PROTOCOL	HY363	12/08/2025	3.999	3.999			0.042	0.033	NA	NA	0.157	
788	TeSeE SAP COMBI	4H0095	04/07/2025	3.408	9.000			0.023	0.017	0.020	0.022	0.230	
829	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	AY185	08/01/2025	4.051	4.020			0.023	0.021	N/A	N/A	0.142	
862	TeSeE SAP COMBI	4J0094	02/05/2025	2.038	1.983			0.008	0.008	0.013	0.006	0.196-0.217	
892	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	CY816	21/03/2025	3.418	3.544			0.020	0.023	0.023	0.020	0.142	Robot

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 3 di 15	DATA:



TSE EURL
Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta – Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RB24	Report No. 01/2025	Distribution date: 02 nd -03 rd December 2024
--	--------------------	---

904	TeSeE SAP COMBI	4C0093	14/01/2025	3.082	3.013			0.013	0.009	0.015	0.009	0.222	
914	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	EY 700	05/06/2025	3.999	3.999			0.021	0.031	'---	'---	0.146	
922	IDEXX HerdChek BSE/Scrapie Antigen Test-ULTRASHORT PROTOCOL	BY445	28/02/2025	3.660	3.700			0.034	0.034	X	X	0.154	

Table 2. General information about '1st alternative' methods applied by the participants in the RB24 PT

NRL ID code	Test Method used	Kit Batch Number	Kit Expiry Date	Positive Control 1 result	Positive Control 2 result	Positive Control 3 result	Positive Control 4 result	Negative Control 1 result	Negative Control 2 result	Negative Control 3 result	Negative Control 4 result	Positive Cut Off value	Methods details
132	TeSeE SAP COMBI	4C0093	14/01/2025	3.087	2.333			0.018	0.015	0.015	0.016	0.226	
135	TeSeE SAP COMBI	4C0093	14/01/2025	2.760	2.732			0.010	0.012	0.012	0.011	0.221	
740	TeSeE SAP COMBI	4J0094	02/05/2025	1.692	1.715			0.008	0.009	0.009	0.009	0.218	
904	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	CY816	21/03/2025	3.5/Over	3.442			0.029	0.022	No negative control was used	3rd negative control was used	0.145	

Table 3. General information about '2nd alternative' methods applied by the participants in the RB24 PT

NRL ID code	Test Method used	Kit Batch Number	Kit Expiry Date	Positive Control1 result	Positive Control2 result	Positive Control3 result	Positive Control4 result	Negative Control1 result	Negative Control2 result	Negative Control3 result	Negative Control4 result	Positive Cut Off value	Methods details
740	BetaPrion BSE EIA Test Kit	0624-01	30/06/2025	2.679	2.679			0.012	0.012	0.012	0.010	0.200	
904	IDEXX HerdChek BSE/Scrapie Antigen Test-ULTRASHORT PROTOCOL	CY816	21/03/2025	3.5/Over	3.5/Over			0.022	0.013	No negative control was used	3rd negative control was used	0.137	

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 4 di 15	DATA:



TSE EURL
Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta - Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RB24	Report No. 01/2025	Distribution date: 02 nd -03 rd December 2024
--	--------------------	---

Table 4. Results submitted by the participants in the RB24 PT applying 'routine' methods

NRL ID code	Test Method used	Run	RB2401 Interpr.	RB2401 OD results	RB2402 Interpr.	RB2402 OD results	RB2403 Interpr.	RB2403 OD results	RB2404 Interpr.	RB2404 OD results	RB2405 Interpr.	RB2405 OD results	RB2406 Interpr.	RB2406 OD results	RB2407 Interpr.	RB2407 OD results
Intended --->			Negative Bovine Tissue Homogenate		Classical Positive Bovine Tissue Homogenate - Dilution 1:4		Negative Bovine Tissue Homogenate		Negative Bovine Tissue Homogenate		Negative Bovine Tissue Homogenate		Negative Bovine Tissue Homogenate		Negative Bovine Tissue Homogenate	
022	TeSeE SAP COMBI	1st	negative	0.010	Medium positive	0.818	negative	0.011	negative	0.010	negative	0.010	negative	0.011	negative	0.011
022	TeSeE SAP COMBI	2nd	negative	0.011	Medium positive	0.774	negative	0.011	negative	0.010	negative	0.010	negative	0.011	negative	0.011
027	TeSeE SAP COMBI	1st	negative	0.022	Strong positive	2.515	negative	0.014	negative	0.020	negative	0.020	negative	0.020	negative	0.020
027	TeSeE SAP COMBI	2nd	negative	0.019	Strong positive	2.472	negative	0.028	negative	0.007	negative	0.013	negative	0.017	negative	0.019
091	TeSeE SAP COMBI	1st	negative	0.018	Medium positive	1.506	negative	0.018	negative	0.026	negative	0.021	negative	0.019	negative	0.017
091	TeSeE SAP COMBI	2nd	negative	0.020	Medium positive	1.718	negative	0.017	negative	0.019	negative	0.026	negative	0.020	negative	0.021
132	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	1st	negative	0.032	Strong positive	3.579	negative	0.026	negative	0.040	negative	0.024	negative	0.027	negative	0.031
132	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	2nd	negative	0.032	Strong positive	3.280	negative	0.021	negative	0.025	negative	0.023	negative	0.029	negative	0.041

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 5 di 15	DATA:



TSE EURL

Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta – Turin

[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RB24 Report No. 01/2025 Distribution date: 02nd-03rd December 2024

135	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHOR T PROTOCOL	1st	negative	0.033	Medium positive	3.349	negative	0.031	negative	0.031	negative	0.034	negative	0.043	negative	0.031
135	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHOR T PROTOCOL	2nd	negative		negative		negative		negative		negative		negative		negative	
181	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHOR T PROTOCOL	1st	negative	0.024	Strong positive	3.999	negative	0.026	negative	0.024	negative	0.027	negative	0.022	negative	0.022
181	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHOR T PROTOCOL	2nd	negative	0.024	Strong positive	3.999	negative	0.024	negative	0.020	negative	0.027	negative	0.028	negative	0.022
192	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHOR T PROTOCOL	1st	negative	0.023	Strong positive	3.557	negative	0.021	negative	0.021	negative	0.021	negative	0.022	negative	0.022

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 6 di 15	DATA:



TSE EURL

Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta - Turin

[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RB24 Report No. 01/2025 Distribution date: 02nd-03rd December 2024

192	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHOR T PROTOCOL	2nd	negative	0.024	Strong positive	3.420	negative	0.023	negative	0.021	negative	0.020	negative	0.023	negative	0.021
291	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	1st	negative	0.023	Strong positive	2.726	negative	0.019	negative	0.018	negative	0.019	negative	0.018	negative	0.019
291	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	2nd	negative	0.025	Strong positive	2.643	negative	0.017	negative	0.020	negative	0.018	negative	0.018	negative	0.019
326	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHOR T PROTOCOL	1st	negative	0.018	Strong positive	2.777	negative	0.016	negative	0.015	negative	0.020	negative	0.007	negative	0.008
326	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHOR T PROTOCOL	2nd	negative	0.029	Strong positive	2.892	negative	0.019	negative	0.017	negative	0.027	negative	0.041	negative	0.014
329	TeSeE SAP COMBI	1st	negative	0.017	Medium positive	0.659	negative	0.012	negative	0.013	negative	0.017	negative	0.012	negative	0.014

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 7 di 15	DATA:



TSE EURL
Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta - Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RB24	Report No. 01/2025	Distribution date: 02 nd -03 rd December 2024
--	--------------------	---

329	TeSeE SAP COMBI	2nd	negative	0.008	Medium positive	0.390	negative	0.013	negative	0.015	negative	0.012	negative	0.009	negative	0.010
358	TeSeE SAP COMBI	1st	negative	0.018	Weak positive	0.660	negative	0.022	negative	0.020	negative	0.022	negative	0.021	negative	0.018
358	TeSeE SAP COMBI	2nd	negative	0.019	Weak positive	0.755	negative	0.020	negative	0.016	negative	0.017	negative	0.018	negative	0.014
390 433 (*)	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	1st	negative	0.036	Strong positive	2.905	negative	0.028	negative	0.027	negative	0.030	negative	0.030	negative	0.028
390 433 (*)	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	2nd	negative	0.026	Strong positive	2.866	negative	0.025	negative	0.024	negative	0.027	negative	0.027	negative	0.026
051 397 (*)	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHORT PROTOCOL	1st	negative	0.018	Strong positive	3.523	negative	0.020	negative	0.009	negative	0.020	negative	0.019	negative	0.021
051 397 (*)	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHORT PROTOCOL	2nd	negative	0.020	Strong positive	3.535	negative	0.023	negative	0.024	negative	0.018	negative	0.017	negative	0.019
415	IDEXX HerdChek BSE/Scrapie	1st	negative	0.016	Strong positive	3.709	negative	0.022	negative	0.016	negative	0.019	negative	0.017	negative	0.020

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 8 di 15	DATA:



TSE EURL
Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta - Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RB24	Report No. 01/2025	Distribution date: 02 nd -03 rd December 2024
--	--------------------	---

	Antigen Test- SHORT PROTOCOL															
415	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	2nd	negative	0.019	Strong positive	3.137	negative	0.020	negative	0.022	negative	0.017	negative	0.017	negative	0.017
595	TeSeE SAP COMBI	1st	negative	0.018	Strong positive	1.769	negative	0.018	negative	0.018	negative	0.017	negative	0.012	negative	0.015
595	TeSeE SAP COMBI	2nd	negative	0.014	Strong positive	2.279	negative	0.016	negative	0.045	negative	0.015	negative	0.014	negative	0.016
596	TeSeE SAP COMBI	1st	negative	0.011	Medium positive	1.385	negative	0.012	negative	0.013	negative	0.014	negative	0.013	negative	0.014
596	TeSeE SAP COMBI	2nd	negative	0.013	Medium positive	1.207	negative	0.012	negative	0.012	negative	0.013	negative	0.014	negative	0.014
651	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	1st	negative	0.047	Medium positive	2.001	negative	0.054	negative	0.031	negative	0.035	negative	0.032	negative	0.036
651	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	2nd	negative	0.125	Strong positive	2.930	negative	0.034	negative	0.033	negative	0.033	negative	0.062	negative	0.039
672	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	1st	negative	0.019	Strong positive	3.022	negative	0.019	negative	0.018	negative	0.024	negative	0.024	negative	0.022

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 9 di 15	DATA:



TSE EURL

Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta - Turin

[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RB24	Report No. 01/2025	Distribution date: 02 nd -03 rd December 2024
--	--------------------	---

672	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	2nd	negative	0.019	Strong positive	3.159	negative	0.018	negative	0.018	negative	0.024	negative	0.020	negative	0.020
703	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHOR T PROTOCOL	1st	negative	0.003	Strong positive	3.705	negative	0.013	negative	0.016	negative	0.005	negative	0.005	negative	0.015
703	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHOR T PROTOCOL	2nd	negative	0.001	Strong positive	overflow	negative	0.005	negative	0.011	negative	0.003	negative	0.007	negative	0.005
740	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHOR T PROTOCOL	1st	negative	0.034	Strong positive	3.999	negative	0.036	negative	0.034	negative	0.032	negative	0.039	negative	0.037
740	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHOR T PROTOCOL	2nd	negative	0.033	Strong positive	3.999	negative	0.033	negative	0.035	negative	0.030	negative	0.039	negative	0.034

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 10 di 15	DATA:



TSE EURL
Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta - Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RB24	Report No. 01/2025	Distribution date: 02 nd -03 rd December 2024
--	--------------------	---

788	TeSeE SAP COMBI	1st	negative	0.039	Weak positive	0.358	negative	0.040	negative	0.048	negative	0.040	negative	0.039	negative	0.042
788	TeSeE SAP COMBI	2nd	negative	0.040	Weak positive	0.305	negative	0.038	negative	0.038	negative	0.032	negative	0.039	negative	0.041
829	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	1st	negative	0.023	Strong positive	3.724	negative	0.037	negative	0.021	negative	0.024	negative	0.024	negative	0.025
829	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	2nd	negative	0.025	Strong positive	3.756	negative	0.025	negative	0.021	negative	0.024	negative	0.021	negative	0.023
862	TeSeE SAP COMBI	1st	negative	0.015	Medium positive	2.030	negative	0.017	negative	0.015	negative	0.015	negative	0.017	negative	0.016
862	TeSeE SAP COMBI	2nd	negative	0.016	Medium positive	2.203	negative	0.017	negative	0.015	negative	0.016	negative	0.016	negative	0.018
892	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	1st	negative	0.022	Strong positive	3.419	negative	0.023	negative	0.020	negative	0.031	negative	0.027	negative	0.024
892	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	2nd	negative	0.021	Strong positive	1.955	negative	0.029	negative	0.023	negative	0.032	negative	0.029	negative	0.022
904	TeSeE SAP COMBI	1st	negative	0.012	Medium positive	1.689	negative	0.018	negative	0.012	negative	0.003	negative	0.016	negative	0.010

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 11 di 15	DATA:



TSE EURL
Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta - Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RB24	Report No. 01/2025	Distribution date: 02 nd -03 rd December 2024
--	--------------------	---

904	TeSeE SAP COMBI	2nd	negative	0.014	Medium positive	2.074	negative	0.017	negative	0.015	negative	0.007	negative	0.012	negative	0.014
914	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	1st	negative	0.024	Strong positive	2.754	negative	0.024	negative	0.026	negative	0.024	negative	0.021	negative	0.021
914	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	2nd	negative	0.021	Strong positive	3.496	negative	0.035	negative	0.021	negative	0.024	negative	0.028	negative	0.018
922	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHO RT PROTOCOL	1st	negative	0.032	Strong positive	3.850	negative	0.030	negative	0.036	negative	0.031	negative	0.030	negative	0.032
922	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHO RT PROTOCOL	2nd	negative	0.027	Strong positive	3.764	negative	0.029	negative	0.034	negative	0.030	negative	0.030	negative	0.030

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 12 di 15	DATA:



TSE EURL
Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta - Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RB24	Report No. 01/2025	Distribution date: 02 nd -03 rd December 2024
--	--------------------	---

Table 5. Results submitted by the participants in the RB24 PT applying '1st alternative' methods

NRL ID code	Test Method used	Run	RB2401 Interpr.	RB2401 OD results	RB2402 Interpr.	RB2402 OD results	RB2403 Interpr.	RB2403 OD results	RB2404 Interpr.	RB2404 OD results	RB2405 Interpr.	RB2405 OD results	RB2406 Interpr.	RB2406 OD results	RB2407 Interpr.	RB2407 OD results
Intended --->			Negative Bovine Tissue Homogenate		Classical Positive Bovine Tissue Homogenate - Dilution 1:4		Negative Bovine Tissue Homogenate		Negative Bovine Tissue Homogenate		Negative Bovine Tissue Homogenate		Negative Bovine Tissue Homogenate		Negative Bovine Tissue Homogenate	
132	TeSeE SAP COMBI	1st	negative	0.020	Medium positive	1.398	negative	0.019	negative	0.017	negative	0.020	negative	0.020	negative	0.017
132	TeSeE SAP COMBI	2nd	negative	0.018	Medium positive	0.992	negative	0.018	negative	0.015	negative	0.018	negative	0.016	negative	0.018
135	TeSeE SAP COMBI	1st	negative	0.013	Weak positive	0.225	negative	0.013	negative	0.013	negative	0.014	negative	0.014	negative	0.015
135	TeSeE SAP COMBI	2nd	////	////	////	////	////	////	////	////	////	////	////	////	////	////
740	TeSeE SAP COMBI	1st	negative	0.013	Strong positive	1.559	negative	0.013	negative	0.015	negative	0.009	negative	0.011	negative	0.013
740	TeSeE SAP COMBI	2nd	negative	0.015	Strong positive	1.508	negative	0.014	negative	0.015	negative	0.012	negative	0.012	negative	0.014
904	IDEXX HerdChek BSE/Scrapie Antigen Test-SHORT PROTOCOL	1st	negative	0.027	Strong positive	3.323	negative	0.032	negative	0.025	negative	0.031	negative	0.031	negative	0.057
904	IDEXX HerdChek BSE/Scrapie Antigen Test-SHORT PROTOCOL	2nd	negative	0.026	Strong positive	3.5/over	negative	0.031	negative	0.023	negative	0.022	negative	0.027	negative	0.027

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 13 di 15	DATA:



TSE EURL
Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta - Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RB24 Report No. 01/2025 Distribution date: 02nd-03rd December 2024

Table 6. Results submitted by the participants in the RB24 PT applying '2nd alternative' methods

NRL ID code	Test Method used	Run	RB2401 Interpr.	RB2401 OD results	RB2402 Interpr.	RB2402 OD results	RB2403 Interpr.	RB2403 OD results	RB2404 Interpr.	RB2404 OD results	RB2405 Interpr.	RB2405 OD results	RB2406 Interpr.	RB2406 OD results	RB2407 Interpr.	RB2407 OD results
Intended --->			Negative Bovine Tissue Homogenate		Classical Positive Bovine Tissue Homogenate - Dilution 1:4		Negative Bovine Tissue Homogenate		Negative Bovine Tissue Homogenate		Negative Bovine Tissue Homogenate		Negative Bovine Tissue Homogenate		Negative Bovine Tissue Homogenate	
740	BetaPrion BSE EIA Test Kit	1st	negative	0.016	Strong positive	2.659	negative	0.015	negative	0.021	negative	0.013	negative	0.014	negative	0.032
740	BetaPrion BSE EIA Test Kit	2nd	negative	0.015	Strong positive	2.658	negative	0.016	negative	0.037	negative	0.011	negative	0.025	negative	0.029
904	IDEXX HerdChek BSE/Scrapie Antigen Test-ULTRASHORT PROTOCOL	1st	negative	0.023	Strong positive	3.404	negative	0.025	negative	0.022	negative	0.013	negative	0.023	negative	0.023
904	IDEXX HerdChek BSE/Scrapie Antigen Test-ULTRASHORT PROTOCOL	2nd	negative	0.018	Strong positive	3.384	negative	0.020	negative	0.017	negative	0.013	negative	0.021	negative	0.018

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 14 di 15	DATA:



TSE EURL

Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta – Turin

[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RB24	Report No. 01/2025	Distribution date: 02 nd -03 rd December 2024
--	--------------------	---

Comments provided by Maria Mazza, IZSPLV - Turin (maria.mazza@izsplv.it)

RB2401 – Negative bovine tissue homogenate

All laboratories reported the correct results.

RB2402 – Classical positive bovine tissue homogenate - Dilution 1:4

All laboratories reported the correct results.

RB2403 - Negative bovine tissue homogenate

All laboratories reported the correct results.

RB2404 - Negative bovine tissue homogenate

All laboratories reported the correct results.

RB2405 - Negative bovine tissue homogenate

All laboratories reported the correct results.

RB2406 - Negative bovine tissue homogenate

All laboratories reported the correct results.

RB2407 - Negative bovine tissue homogenate

All laboratories reported the correct results.

✚ (*) The laboratories 051 and 390 entered incorrect ID codes instead of the correct codes 397 and 433, respectively.

✚ The laboratory 135 did not analyse the samples in duplicate, as requested.

✚ The laboratories 415 and 390 (433) reported optical densities results of the samples in duplicate but in two different analysis sessions and not in the same session. It is recommended for the upcoming PTs to analyse the samples in duplicate in the same run to assess the repeatability.

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 15 di 15	DATA:



TSE EURL
Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta – Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RB24	Report No. 01/2025	Distribution date: 02 nd -03 rd December 2024
--	--------------------	---

Conclusions provided by Maria Mazza, IZSPLV - Turin (maria.mazza@izsplv.it)

In summary all laboratories have passed this PT round. Well done.

Please remind that, according to Reg. 999/2001 (EC), rapid tests may be used for both primary screening of suspect cases and, if inconclusive or positive, for subsequent confirmation ...*omissis*.... provided that:
...*omissis*....

The information contained in this Report must be considered confidential

Date: 24th February 2025

Dr. Elena Bozzetta, Coordinator of TSE EURL

Dr. Giuseppe Ru, Director of TSE EURL

END OF REPORT

Copia controllata