

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 1 di 16	DATA:



TSE EURL

Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta – Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RS24	Report No. 02/2025	Distribution date: 2 nd -3 rd December 2024
--	--------------------	---

Table 1. General information about 'routine' methods applied by the participants in the RS24 PT

NRL ID code	Test Method used	Kit Batch Number	Kit Expiry Date	Positive Control1 result	Positive Control2 result	Positive Control3 result	Positive Control4 result	Negative Control1 result	Negative Control2 result	Negative Control3 result	Negative Control4 result	Positive Cut Off value	Methods details
105	IDEXX HerdChek BSE/Scrapie Antigen Test-SHORT PROTOCOL	BY445	28/02/2025	3.873	4.076			0.031	0.029	'-	'-	0.210	
107	IDEXX HerdChek BSE/Scrapie Antigen Test-SHORT PROTOCOL	AY185	08/01/2025	3.350	3.358			0.033	0.033	0.033	0.033	0.213	
193	TeSeE SAP COMBI	4H0095	04/07/2025	3.249	3.268			0.007	0.011	0.008	0.008	Negative cut-off < 0.197; Positive cut-off >= 0.219	
200	IDEXX HerdChek BSE/Scrapie Antigen Test-ULTRASHORT PROTOCOL	BY445	28/02/2025	1.878	2.000			0.039	0.041	'-	'-	0.220	
234	TeSeE SAP COMBI	4H0095	04/07/2025	* **	* **			0.008	0.009	0.008	0.009	0.219	All samples were tested twice.
244	IDEXX HerdChek BSE/Scrapie Antigen Test-SHORT PROTOCOL	EY700	05/06/2025	3.999	3.999			0.033	0.033	0.031	0.031	0.212	
254	TeSeE SAP COMBI	4C0093	14/01/2025	3.649	3.633			0.010	0.010	0.011	0.006	0.197	
257	IDEXX HerdChek BSE/Scrapie Antigen Test-SHORT PROTOCOL	CY816 and EY700	21/03/2025	3.396	3.540			0.030	0.032	0.061	0.064	0,211 and 0,243	We used 2 different kit batches, ran 05.12.2024 CY816, cut off 0,211, and ran 19.12.2024 cut off 0,243. The results are in 2

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 2 di 16	DATA:



TSE EURL

Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta – Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RS24	Report No. 02/2025	Distribution date: 2 nd -3 rd December 2024
--	--------------------	---

													separate files, I sent the first run in the RB24 file, and send the second run in RS24, sorry for that.
275	IDEXX HerdChek BSE/Scrapie Antigen Test-SHORT PROTOCOL	EY700	05/06/2025	3.653	4.008			0.026	0.025	0.035	0.031	0.213	
278	TeSeE SAP COMBI	4J0094	02/05/2025	3.557	3.530			0.014	0.015	0.016	0.011	0.224	
338	TeSeE SAP COMBI	4C0093	14/01/2025	3.096	2.976			0.013	0.008	0.016	0.010	0.222	
389	IDEXX HerdChek BSE/Scrapie Antigen Test-SHORT PROTOCOL	EY700	05/06/2025	3.915	3.887			0.043	0.034	0.043	0.034	0.218	
447	IDEXX HerdChek BSE/Scrapie Antigen Test-ULTRASHORT PROTOCOL	BY445	28/02/2025	3.857	3.746			0.035	0.035	X	X	0.215	
467	TeSeE SAP COMBI	4H0095	04/07/2025	3.408	9.000			0.023	0.017	0.020	0.022	0.230	
468	TeSeE SAP COMBI	4C0093	14/01/2025	3.207	3.215			0.004	0.007	0.007	0.006	0.216	
051 510 (*)	IDEXX HerdChek BSE/Scrapie Antigen Test-ULTRASHORT PROTOCOL	EY700	05/06/2025	3.589	3.568			0.023	0.023			0.203	
515	IDEXX HerdChek BSE/Scrapie Antigen Test-SHORT PROTOCOL	AY185	08/01/2025	3.389	3.430			0.030	0.029			0.209	
543	IDEXX HerdChek BSE/Scrapie Antigen Test-SHORT PROTOCOL	AY185	08/01/2025	4.011	4.011			0.032	0.035	N/A	N/A	0.214	

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 3 di 16	DATA:



TSE EURL

Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta – Turin

Istituto Superiore di Sanità - Rome

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RS24	Report No. 02/2025	Distribution date: 2 nd -3 rd December 2024
--	--------------------	---

613	TeSeE SAP COMBI	4J0094	02/05/2025	2.080	2.192			0.008	0.009	0.009	0.008	0.218	
639	TeSeE SAP COMBI	4J0094	02/05/2025	1.721	1.688			0.006	0.006	0.008	0.007	0.217	
660	TeSeE SAP COMBI	4C0093	14/01/2025	3.072	3.196			0.006	0.006	0.007	0.009	0.217	
718	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	EY 700	05/06/2025	3.999	3.999			0.030	0.029	'---	'---	0.210	
778	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHORT PROTOCOL	HY363	12/08/2025	3.999	3.999			0.056	0.039	NA	NA	0.228	
807	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHORT PROTOCOL	CY816	21/03/2025	>4 (overflow)	>4 (overflow)			0.031	0.033	Na	Na	0.180	background subtracted OD reported for samples
810	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	AY185	08/01/2025	3.533	3.744			0.023	0.021	'-	'-	0.202	
831 583 (*)	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	CY816	21/03/2025	2.922 / 2.913	2.946 / 2.926			0.036 / 0.038	0.041 / 0.039	'-	'-	0.217 / 0.220	
856 896 (*)	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHORT PROTOCOL	EY 700	05/06/2025	3.303	3.125			0.026	0.029	Na	Na	0.206	
927	TeSeE SAP COMBI	4J0094	02/05/2025	2.038	1.983			0.008	0.008	0.013	0.006	0.196-0.217	

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 4 di 16	DATA:



TSE EURL

Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta – Turin

Istituto Superiore di Sanità - Rome

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RS24	Report No. 02/2025	Distribution date: 2 nd -3 rd December 2024
--	--------------------	---

Table 2. General information about '1st alternative' methods applied by the participants in the RS24 PT

NRL ID code	Test Method used	Kit Batch Number	Kit Expiry Date	Positive Control 1 result	Positive Control 2 result	Positive Control 3 result	Positive Control 4 result	Negative Control 1 result	Negative Control 2 result	Negative Control 3 result	Negative Control 4 result	Positive Cut Off value	Methods details
105	TeSeE SAP COMBI	4C0093	14/01/2025	3.087	2.333			0.018	0.015	0.015	0.016	0.226	
200	TeSeE SAP COMBI	4C0093	14/01/2025	2.760	2.732			0.010	0.012	0.012	0.011	0.221	
338	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	CY816	21/03/2025	3.440	3.373			0.032	0.031	No 3rd negative control was used	No 4th negative control was used	0.211	
778	TeSeE SAP COMBI	4J0094	02/05/2025	1.692	1.715			0.008	0.009	0.009	0.009	0.218	

Table 3. General information about '2nd alternative' methods applied by the participants in the RS24 PT

NRL ID code	Test Method used	Kit Batch Number	Kit Expiry Date	Positive Control 1 result	Positive Control 2 result	Positive Control 3 result	Positive Control 4 result	Negative Control 1 result	Negative Control 2 result	Negative Control 3 result	Negative Control 4 result	Positive Cut Off value	Methods details
338	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHORT PROTOCOL	CY816	21/03/2025	3.5/Over	3.377			0.032	0.029	No 3rd negative control was used	No 4th negative control was used	0.210	

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 5 di 16	DATA:



TSE EURL

Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta – Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RS24	Report No. 02/2025	Distribution date: 2 nd -3 rd December 2024
--	--------------------	---

Table 4. Results submitted by the participants in the RS24 PT applying 'routine' methods

NRL ID code	Test Method used	Run	RS2401 Interpr.	RS2401 OD results	RS2402 Interpr.	RS2402 OD results	RS2403 Interpr.	RS2403 OD results	RS2404 Interpr.	RS2404 OD results	RS2405 Interpr.	RS2405 OD results	RS2406 Interpr.	RS2406 OD results	RS2407 Interpr.	RS2407 OD results
Intended --->			TSE Negative Ovine Tissue Homogenate		TSE Negative Ovine Tissue Homogenate		TSE Negative Ovine Tissue Homogenate		Classical Scrapie Positive Ovine Tissue Homogenate – Dilution 1:10		TSE Negative Ovine Tissue Homogenate		Classical Scrapie Positive Caprine Tissue Homogenate- Dilution 1:4		TSE Negative ovine Tissue Homogenate	
105	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	1st	negative	0.034	negative	0.033	negative	0.030	medium positive	0.931	negative	0.031	Strong positive	3.559	negative	0.025
105	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	2nd	negative	0.033	negative	0.036	negative	0.028	medium positive	0.917	negative	0.030	Strong positive	3.332	negative	0.023
107	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	1st	negative	0.036	negative	0.152	negative	0.101	medium positive	1.708	negative	0.066	Strong positive	3.142	negative	0.098
107	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	2nd	negative	0.128	negative	0.150	negative	0.068	medium positive	1.114	negative	0.091	Medium positive	2.092	negative	0.066
193	TeSeE SAP COMBI	1st	negative	0.023	negative	0.020	negative	0.020	Weak positive	0.707	negative	0.017	Strong positive	2.495	negative	0.015

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 6 di 16	DATA:



TSE EURL

Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta – Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RS24	Report No. 02/2025	Distribution date: 2 nd -3 rd December 2024
--	--------------------	---

193	TeSeE SAP COMBI	2nd	negative	0.018	negative	0.021	negative	0.019	Weak positive	0.893	negative	0.018	Strong positive	2.569	negative	0.017
200	IDEXX HerdChek BSE/Scrapie Antigen Test-ULTRASHORT PROTOCOL	1st	negative	0.047	negative	0.046	negative	0.044	Weak positive	1.449	negative	0.047	Medium positive	3.332	negative	0.046
200	IDEXX HerdChek BSE/Scrapie Antigen Test-ULTRASHORT PROTOCOL	2nd	negative		negative		negative		negative		negative		negative		negative	
234	TeSeE SAP COMBI	1st	negative	0.021	negative	0.021	negative	0.021	Medium positive	1.015	negative	0.020	Strong positive	* **	negative	0.019
234	TeSeE SAP COMBI	2nd	negative	0.014	negative	0.019	negative	0.021	Medium positive	1.181	negative	0.019	Strong positive	* **	negative	0.021
244	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	1st	negative	0.034	negative	0.031	negative	0.038	Medium positive	2.973	negative	0.030	Strong positive	3.999	negative	0.036
244	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	2nd	negative	0.036	negative	0.035	negative	0.034	Medium positive	2.975	negative	0.035	Strong positive	3.999	negative	0.035
254	TeSeE SAP COMBI	1st	negative	0.010	negative	0.007	negative	0.010	Weak positive	0.388	negative	0.017	Medium positive	0.620	negative	0.142

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 7 di 16	DATA:



TSE EURL
Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta – Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RS24	Report No. 02/2025	Distribution date: 2 nd -3 rd December 2024
--	--------------------	---

254	TeSeE SAP COMBI	2nd	negative	0.011	negative	0.012	negative	0.015	Weak positive	0.291	negative	0.017	Medium positive	1.628	negative	0.015
257	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	1st	negative	0.039	negative	0.044	negative	0.040	Medium positive	1.621	negative	0.050	Strong positive	3.645	negative	0.049
257	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	2nd	negative	0.078	negative	0.072	negative	0.063	Medium positive	1.167	negative	0.067	Strong positive	3.210	negative	0.065
275	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	1st	negative	0.035	negative	0.030	negative	0.041	Medium positive	2.551	negative	0.030	Strong positive	4.100	negative	0.033
275	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	2nd	negative	0.027	negative	0.026	negative	0.030	Medium positive	1.411	negative	0.023	Strong positive	3.572	negative	0.032
278	TeSeE SAP COMBI	1st	negative	0.016	negative	0.016	negative	0.018	Weak positive	0.727	negative	0.017	Weak positive	0.933	negative	0.017
278	TeSeE SAP COMBI	2nd	negative	0.017	negative	0.019	negative	0.019	Weak positive	0.575	negative	0.018	Weak positive	0.910	negative	0.013
338	TeSeE SAP COMBI	1st	negative	0.016	negative	0.012	negative	0.018	Medium positive	1.080	negative	0.005	Medium positive	2.162	negative	0.017
338	TeSeE SAP COMBI	2nd	negative	0.016	negative	0.005	negative	0.015	Medium positive	1.131	negative	0.005	Medium positive	2.273	negative	0.020

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 8 di 16	DATA:



TSE EURL
Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta – Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RS24	Report No. 02/2025	Distribution date: 2 nd -3 rd December 2024
--	--------------------	---

389	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	1st	negative	0.030	negative	0.031	negative	0.032	Medium positive	1.678	negative	0.037	Strong positive	3.842	negative	0.042
389	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	2nd	negative	0.030	negative	0.028	negative	0.034	Medium positive	1.581	negative	0.039	Strong positive	3.746	negative	0.044
447	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHOR T PROTOCOL	1st	negative	0.043	negative	0.040	negative	0.036	Medium positive	1.307	negative	0.035	Strong positive	3.527	negative	0.037
467	TeSeE SAP COMBI	2nd	negative	0.037	negative	0.036	negative	0.036	Medium positive	1.166	negative	0.035	Strong positive	3.797	negative	0.038
467	TeSeE SAP COMBI	2nd	negative	0.039	negative	0.039	negative	0.037	Weak positive	1.004	negative	0.035	Strong positive	2.404	negative	0.040
468	TeSeE SAP COMBI	1st	negative	0.014	negative	0.011	negative	0.015	Medium positive	1.234	negative	0.014	Strong positive	3.456	negative	0.014
468	TeSeE SAP COMBI	2nd	negative	0.013	negative	0.014	negative	0.013	Medium positive	1.221	negative	0.013	Strong positive	3.422	negative	0.012
051 510 (*)	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHOR T PROTOCOL	1st	negative	0.016	negative	0.026	negative	0.026	Medium positive	1.664	negative	0.027	Strong positive	3.523	negative	0.026

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 9 di 16	DATA:



TSE EURL

Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta – Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RS24	Report No. 02/2025	Distribution date: 2 nd -3 rd December 2024
--	--------------------	---

051 510 (*)	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHOR T PROTOCOL	2nd	negative	0.025	negative	0.026	negative	0.028	Medium positive	1.599	negative	0.024	Strong positive	3.512	negative	0.035
515	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	1st	negative	0.061	negative	0.046	negative	0.035	Medium positive	2.309	negative	0.035	Strong positive	3.794	negative	0.038
515	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	2nd	negative	0.049	negative	0.044	negative	0.044	Medium positive	2.454	negative	0.039	Strong positive	3.854	negative	0.043
543	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	1st	negative	0.044	negative	0.051	negative	0.042	Strong positive	2.597	negative	0.033	Strong positive	4.045	negative	0.032
543	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	2nd	negative	0.038	negative	0.041	negative	0.043	Strong positive	2.730	negative	0.036	Strong positive	4.016	negative	0.020
613	TeSeE SAP COMBI	1st	negative	0.013	negative	0.014	negative	0.014	Weak positive	0.301	negative	0.016	Medium positive	1.052	negative	0.014

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 10 di 16	DATA:



TSE EURL
Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta – Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RS24	Report No. 02/2025	Distribution date: 2 nd -3 rd December 2024
--	--------------------	---

613	TeSeE SAP COMBI	2nd	negative	0.017	negative	0.014	negative	0.017	Medium positive	0.736	negative	0.013	Strong positive	2.489	negative	0.014
639	TeSeE SAP COMBI	1st	negative	0.013	negative	0.014	negative	0.012	Weak positive	0.408	negative	0.014	Medium positive	1.491	negative	0.013
639	TeSeE SAP COMBI	2nd	negative	0.013	negative	0.012	negative	0.011	Weak positive	0.448	negative	0.012	Medium positive	1.236	negative	0.014
660	TeSeE SAP COMBI	1st	negative	0.026	negative	0.007	negative	0.010	Weak positive	0.801	negative	0.028	Strong positive	2.570	negative	0.031
660	TeSeE SAP COMBI	2nd	negative	0.025	negative	0.028	negative	0.005	Weak positive	0.850	negative	0.014	Strong positive	2.244	negative	0.026
718	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	1st	negative	0.029	negative	0.037	negative	0.036	Medium positive	1.059	negative	0.034	Strong positive	2.820	negative	0.034
718	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	2nd	negative	0.028	negative	0.037	negative	0.030	Medium positive	2.070	negative	0.033	Strong positive	3.592	negative	0.029
778	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHORT PROTOCOL	1st	negative	0.038	negative	0.043	negative	0.065	Medium positive	1.481	negative	0.035	Strong positive	3.999	negative	0.071
778	IDEXX HerdChek BSE/Scrapie Antigen Test-	2nd	negative	0.041	negative	0.038	negative	0.039	Medium positive	1.456	negative	0.035	Strong positive	3.999	negative	0.043

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 11 di 16	DATA:



TSE EURL

Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta – Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RS24	Report No. 02/2025	Distribution date: 2 nd -3 rd December 2024
--	--------------------	---

	ULTRASHOR T PROTOCOL															
807	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHOR T PROTOCOL	1st	negative	0.004	negative	0.009	negative	0.003	Medium positive	1.695	negative	0.008	Strong positive	overflow	negative	0.009
807	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHOR T PROTOCOL	2nd	negative	0.005	negative	0.021	negative	0.002	Medium positive	1.903	negative	0.010	Strong positive	overflow	negative	0.004
810	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	1st	negative	0.029	negative	0.031	negative	0.030	Medium positive	1.222	negative	0.028	Strong positive	3.795	negative	0.028
810	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	2nd	negative	0.031	negative	0.027	negative	0.030	Medium positive	1.226	negative	0.027	Strong positive	3.785	negative	0.029
831 583 (*)	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	1st	negative	0.046	negative	0.043	negative	0.044	Strong positive	2.540	negative	0.041	Strong positive	2.895	negative	0.045

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 12 di 16	DATA:



TSE EURL

Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta – Turin

[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RS24	Report No. 02/2025	Distribution date: 2 nd -3 rd December 2024
--	--------------------	---

831 583 (*)	IDEXX HerdChek BSE/Scrapie Antigen Test- SHORT PROTOCOL	2nd	negative	0.053	negative	0.045	negative	0.047	Strong positive	2.759	negative	0.053	Strong positive	2.911	negative	0.049
856 896 (*)	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHOR T PROTOCOL	1st	negative	0.030	negative	0.055	negative	0.028	Medium positive	1.543	negative	0.037	Strong positive	3.306	negative	0.031
856 896 (*)	IDEXX HerdChek BSE/Scrapie Antigen Test- ULTRASHOR T PROTOCOL	2nd	negative	0.037	negative	0.076	negative	0.029	Medium positive	1.857	negative	0.042	Strong positive	3.425	negative	0.032
927	TeSeE SAP COMBI	1st	negative	0.016	negative	0.015	negative	0.018	Weak positive	1.110	negative	0.015	Medium positive	2.729	negative	0.018
927	TeSeE SAP COMBI	2nd	negative	0.015	negative	0.017	negative	0.018	Weak positive	1.054	negative	0.015	Medium positive	2.564	negative	0.016

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 13 di 16	DATA:



TSE EURL
Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta – Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RS24	Report No. 02/2025	Distribution date: 2 nd -3 rd December 2024
--	--------------------	---

Table 5. Results submitted by the participants in the RS24 PT applying '1st alternative' methods

NRL ID code	Test Method used	Run	RS2401 Interpr.	RS2401 OD results	RS2402 Interpr.	RS2402 OD results	RS2403 Interpr.	RS2403 OD results	RS2404 Interpr.	RS2404 OD results	RS2405 Interpr.	RS2405 OD results	RS2406 Interpr.	RS2406 OD results	RS2407 Interpr.	RS2407 OD results
Intended ---->			TSE Negative Ovine Tissue Homogenate		TSE Negative Ovine Tissue Homogenate		TSE Negative Ovine Tissue Homogenate		Classical Scrapie Positive Ovine Tissue Homogenate – Dilution 1:10		TSE Negative Ovine Tissue Homogenate		Classical Scrapie Positive Caprine Tissue Homogenate – Dilution 1:4		TSE Negative Ovine Tissue Homogenate	
105	TeSeE SAP COMBI	1st	negative	0.021	negative	0.019	negative	0.018	Medium positive	0.850	negative	0.022	Strong positive	1.825	negative	0.020
105	TeSeE SAP COMBI	2nd	negative	0.020	negative	0.019	negative	0.017	Medium positive	0.664	negative	0.021	Strong positive	2.606	negative	0.020
200	TeSeE SAP COMBI	1st	negative	0.014	negative	0.015	negative	0.013	Weak positive	0.236	negative	0.014	Medium positive	0.605	negative	0.016
200	TeSeE SAP COMBI	2nd	negative		negative		negative		negative		negative		negative		negative	
338	IDEXX HerdChek BSE/Scrapie Antigen Test-SHORT PROTOCOL	1st	negative	0.039	negative	0.033	negative	0.042	Medium positive	2.437	negative	0.038	Strong positive	3.5/Over	negative	0.066
338	IDEXX HerdChek BSE/Scrapie Antigen Test-SHORT PROTOCOL	2nd	negative	0.040	negative	0.034	negative	0.037	Medium positive	2.238	negative	0.036	Strong positive	3.5/Over	negative	0.054
778	TeSeE SAP COMBI	1st	negative	0.012	negative	0.016	negative	0.020	Strong positive	2.333	negative	0.010	Strong positive	2.628	negative	0.013

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 14 di 16	DATA:



TSE EURL

Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta – Turin
[Istituto Superiore di Sanità - Rome](#)

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RS24	Report No. 02/2025	Distribution date: 2 nd -3 rd December 2024
--	--------------------	---

778	TeSeE SAP COMBI	2nd	negative	0.015	negative	0.018	negative	0.021	Strong positive	2.335	negative	0.013	Strong positive	2.613	negative	0.013
-----	-----------------	-----	----------	-------	----------	-------	----------	-------	-----------------	-------	----------	-------	-----------------	-------	----------	-------

Table 6. Results submitted by the participants in the RS24 PT applying '2nd alternative' methods

NRL ID code	Test Method used	Run	RS2401 Interpr.	RS2401 OD results	RS2402 Interpr.	RS2402 OD results	RS2403 Interpr.	RS2403 OD results	RS2404 Interpr.	RS2404 OD results	RS2405 Interpr.	RS2405 OD results	RS2406 Interpr.	RS2406 OD results	RS2407 Interpr.	RS2407 OD results
Intended --->			TSE Negative Ovine Tissue Homogenate		TSE Negative Ovine Tissue Homogenate		TSE Negative Ovine Tissue Homogenate		Classical Scrapie Positive Ovine Tissue Homogenate – Dilution 1:10		TSE Negative Ovine Tissue Homogenate		Classical Scrapie Positive Caprine Tissue Homogenate – Dilution 1:4		TSE Negative Ovine Tissue Homogenate	
338	IDEXX HerdChek BSE/Scrapie Antigen Test-ULTRASHO RT PROTOCOL	1st	negative	0.033	negative	0.026	negative	0.037	Medium positive	2.484	negative	0.032	Strong positive	3.431	negative	0.055
338	IDEXX HerdChek BSE/Scrapie Antigen Test-ULTRASHO RT PROTOCOL	2nd	negative	0.028	negative	0.023	negative	0.032	Strong positive	2.638	negative	0.026	Strong positive	3.370	negative	0.059

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 15 di 16	DATA:



TSE EURL

Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta – Turin

[Istituto Superiore di Sanità](#) - Rome

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RS24	Report No. 02/2025	Distribution date: 2 nd -3 rd December 2024
--	--------------------	---

Comments provided by Maria Mazza, IZSPLVA - Turin (maria.mazza@izsplv.it)

RS2401 - TSE Negative ovine tissue homogenate

All laboratories reported the correct result.

RS2402 - TSE Negative ovine tissue homogenate

All laboratories reported the correct result.

RS2403 - TSE Negative ovine tissue homogenate

All laboratories reported the correct result.

RS2404 - Classical Scrapie positive ovine tissue homogenate - Dilution 1:10

All laboratories reported the correct result.

RS2405 - TSE Negative ovine tissue homogenate

All laboratories reported the correct result.

RS2406 - Classical Scrapie positive caprine tissue homogenate - Dilution 1:4

All laboratories reported the correct result.

RS2407 - TSE Negative ovine tissue homogenate

All laboratories reported the correct result.

✚ (*) The laboratories 051, 856 and 831 entered incorrect ID codes instead of the correct codes 510, 896 and 583, respectively.

✚ The laboratory 856 (896) did not report the batch number of the kit used.

✚ The laboratory 200 did not analyse the samples in duplicate, as requested.

✚ The laboratories 244, 257, 275, 831 (583) and 927 reported optical densities results of the samples in duplicate but in two different analysis sessions and not in the same session. The laboratory 257 also used two different batches for the analysis.

It is recommended for the upcoming PTs to analyse the samples in duplicate in the same run to assess the repeatability.

Results Tabulation for Screening Rapid test

MOD. N°	40IZ476	ED/REV. 00/01	STRUTTURA/UO:
PGS. N°	10DG010	PAGINA 16 di 16	DATA:



TSE EURL

Istituto Zooprofilattico Sperimentale del
Piemonte, Liguria e Valle d'Aosta – Turin

[Istituto Superiore di Sanità](#) - Rome

PROFICIENCY TESTING FOR VETERINARY LABORATORIES

Results Tabulation for Screening Rapid test - Scheme: RS24	Report No. 02/2025	Distribution date: 2 nd -3 rd December 2024
--	--------------------	---

Conclusions provided by Maria Mazza, IZSPLVA - Turin (maria.mazza@izsplv.it)

In summary all laboratories have passed this PT round. Well done.

Please remind that, according to Reg. 999/2001 (EC), rapid tests may be used for both primary screening of suspect cases and, if inconclusive or positive, for subsequent confirmation ...*omissis*.... provided that:

...*omissis*....

The information contained in this Report must be considered confidential

Date: 24th February 2025

Dr. Elena Bozzetta, Coordinator of TSE EURL

Dr. Giuseppe Ru, Director of TSE EURL

END OF REPORT

Copia controllata